



Certificate of Conformity

VWR CARDBOARD CRYOBBOXES 100 mm HIGH

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	WxDxH ext. (mm)	Colour	Height (mm)	Pk
479-1392	136x136x100	White	100	1
479-1393	136x136x100	Blue	100	1
479-1394	136x136x100	Red	100	1
479-1395	136x136x100	Green	100	1
479-1396	136x136x100	Yellow	100	1


Material: Cardboard

Other Characteristics:

- For storage of cryovials at low temperatures.
- Waterproof coating.
- Various combinations possible with different dividers (to be ordered separately).
- Available in different heights and 5 colours.
- Freezable to -140 °C



Certificate of Conformity



Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **January 2017**