



Certificate of Conformity

VWR CELL ATTACHMENT TISSUE CULTURE DISHES, TREATED

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No.	Description	Growth Area	Qty/Pk	Qty/Case
734-2814	3.5cm, Treated for Increased Cell Attachment, Sterilized	8.5 cm ²	40	240
734-2815	6.0cm, Treated for Increased Cell Attachment, Sterilized	21.2 cm ²	40	240
734-2816	9.0cm, Treated for Increased Cell Attachment, Sterilized	55 cm ²	40	240
734-2817	10.0cm, Treated for Increased Cell Attachment, Sterilized	60.8 cm ²	40	240
734-2818	15.0cm, Treated for Increased Cell Attachment, Sterilized	143 cm ²	20	80

Materials:


Dish: GPPS (General polystyrene) Colour: Clear

Lid: GPPS(General polystyrene) Colour: Clear



Certificate of Conformity


734-2814





Certificate of Conformity


734-2815





Certificate of Conformity


734-2816





Certificate of Conformity


734-2817





Certificate of Conformity

734-2818



Certificate of Conformity

Other Characteristics:

- Temperature: -20°C~+50°C
- Sterilized by gamma irradiation SAL10⁻⁶(ISO11137)
- Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production.
- Product is latex free.



Certificate of Conformity

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **March 2017**