



## Certificate of Conformity

# VWR CRYOGENIC VIALS WITH INTERNAL THREAD AND SILICONE O-RING

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

*VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document*

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

| Catalog No. | Capacity (ml) | Base         | Ø ext.×H (mm) | Pk    |
|-------------|---------------|--------------|---------------|-------|
| 479-1261    | 1.2           | Freestanding | 12.5×41       | 1.000 |
| 479-1262    | 2.0           | Freestanding | 12.5×49       | 1.000 |
| 479-1263    | 2.0           | Round bottom | 12.5×48       | 1.000 |
| 479-1264    | 4.0           | Freestanding | 12.5×72       | 1.000 |
| 479-1265    | 4.0           | Round bottom | 12.5×70       | 1.000 |
| 479-1266    | 5.0           | Round bottom | 12.5×91       | 1.000 |
| 479-1773    | 5.0           | Freestanding | 12.5×91       | 1.000 |

Material: **PP vial and closure, sterile (SAL 10<sup>-6</sup>)**

Temperature range: **-196°C to +121°C Autoclavable at 121°C, for up to 30 minutes.**

Product compliance:

**Resin conforms to FDA 21 CFR 177.1520**

**Resin conforms to USP Class VI**

**Plastics and colorants are in conformity with CONEG standards**

Other Characteristics:

**Certified RNase, DNase, Pyrogen and DNA free**

**Gamma radiation sterilized.**

**Printed graduations and writing area**

**Clear , no presence of contamination in plastic.**

**Closure verification**

**Leak proof testing in vacuum of 71.3cm/Hg.**

**Gas phase of Liquid nitrogen resistance testing.**

**Disclaimer:** VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **March 2017**