



Certificate of Conformity

VWR GRID DIVIDERS FOR CRYOBBOXES, 100/130 mm

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	Compartments (mm)	Height (mm)	Pk
479-1443	3 x 3	65	1
479-1444	4 x 4	65	1
479-1445	5 x 5	65	1
479-1446	6 x 6	65	1
479-1447	7 x 7	65	1
479-1448	8 x 8	65	1
479-1449	9 x 9	65	1
479-1450	10 x 10	65	1

Material: Cardboard


Other Characteristics:

- Resistance to temperature -140 °C.
- Waterproof coating.
- For storage of cryovials at low temperatures.




Certificate of Conformity

479-1443




479-1444






Certificate of Conformity

479-1445




479-1446






Certificate of Conformity

479-1447




479-1448






Certificate of Conformity

479-1449



479-1450



Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.



Certificate of Conformity

Date of issue: **May 2017**