



Certificate of Conformity

VWR PIPETTE KIT WITH BLUE BULB FOR DROPPING BOTTLES

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No.	Capacity (ml)	Length (cm)	Pk
215-3823	5	5	100
215-3824	10 and 15	6	100
215-3825	20	7	100
215-3826	30	7.5	100
215-3827	50	9	100
215-3828	100	11	100

Material: Soda lime glass, Clear

Other Properties:

- Equipped with a blue rubber bulb and a screw cap made from white PE.
- Good chemical and physical properties.
- Delivered in a plastic bag.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **April 2017**