



Certificate of Conformity

VWR PRECISION BAROMETER, TRACEABLE®

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat No: **628-0041**

Accuracy: **±0.5 %**

Range: **954 - 1073 mbar; 715 - 805 mm Hg**

Resolution: **0.5 mbar; 0.5 mm Hg**

Other Characteristics:

- This barometer with triple scale provides a sensitive measure of atmospheric pressure.
- Corrosion-proof, nickel-chrome.
- Three distinct scales.
- Adjustable pointer to indicate changes.
- Response time to change from one to four hours.
- Compliant with ISO/IEC 17025, ANSI/NCSL Z540-1 and 9001.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **July 2017**