



Certificate of Conformity

VWR SCREW CAPS FOR NARROW MOUTH BOTTLES WITH THREAD

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

VWR Catalog Number	Thread (mm)	For bottles	Pk
215-5643	18	50/100 ml	10
215-5644	25	250/500 ml	10
215-5645	28	1000/2000 ml	10

Material: LDPE, translucent

Other Characteristics:

- Excellent for sample collection, storage and transporting liquids.
- Leakproof, robust and versatile
- Squeezable, containers can be used as dropper or wash bottles
- Raised graduations

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **September 2017**