



Certificate of Conformity

VWR CLEANING BRUSH FOR SMALL VESSELS

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No. : 149-0203

Material : White nylon bristles with wire handle

Head Ø (mm) : 5 - 13

Trimming length (mm) : 75

Total length (mm) : 180

Pk: 5

Other Characteristics:

- These brushes are specially designed for cleaning small tubes, pipes, and a wide range of small containers.
- Rounded tip for gentle but thorough cleaning
- Wire handle for secure grip
- With hanging loop

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **September 2015**