



Certificate of Conformity

VWR DIGITAL STOPWATCH

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Cat.Number : 609-0139

Accuracy : 0.01 sec

Colour : Black

WxH : 63 x 80 mm

Pk: 1

Other Characteristics:

This electronic stopwatch includes a calendar and a clock.

1/100 second chronograph

Alarm function

Memory function (stores one lap time)

CE marked and RoHS complaint

Supplied with one CR2032 battery

Certifications: Tested and calibrated in a triple accredited/registered ISO/IEC 17025:2005, A2LA accredited laboratory (Cert #2448.01), against equipment whose calibration is traceable to NIST.

.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **December 2017**