


THERMOMETER CALIBRATIONS

Can you trust your thermometer?

- In-house calibration facility in Bristol
- UKAS calibration of total and partial immersion liquid in glass thermometers between -35 °c and 650 °c
- UKAS calibration of digital thermometers (system calibrations) between -35 °c and 650 °c
- Best measurement capability of 0.024 °c using our high specification calibration equipment



In-house 3 point UKAS thermometer calibration



In-house 3 point UKAS thermometer calibration. **PRICE:** £80.00 plus carriage (ex VAT). Extra points can be calibrated at £12.00 per extra point if required. You will receive a UKAS approved calibration certificate together with the raw data, providing evidence of the temperature recorded.

Email: lab3service@uk.vwr.com